**Technical Features :**

1. Product size (L x W x H): 219x 64 x 18cm
2. Packing size: 223 x 68x 21cm
3. N.W.: 18kg G.W.: 20kg
4. Load bearing: less than 270kg

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**STRETCHER TC102C: UNIVERSAL RESCUE STRETCHER Use instructions:** -GENERAL USE INFORMATION: This manual must be translated (according to regulations) by the dealer, in the language of the country where the equipment is used. This manual must be read and understood by the user before using the PPE or the rescue device. These instructions for use should be retained throughout the life of the product. Read the instruction manual before use. Contact us for further information; The serial number as indicated below can be found on each appliance and must not be removed or covered. In order to facilitate assistance please indicate or communicate the serial number (SN) or lot number (LOT) on the label. •**WARNINGS:** Follow step by step the instructions and recommendations of the manufacturer in order to guarantee a suitable installation. If in doubt, contact the manufacturer. For your safety, carefully follow the instructions for use, verification, maintenance and storage, described in this leaflet. In case of doubt on the status of the device (traces of oxidation) or after a fall (deformation), it should not be reused and / or be returned to the manufacturer or a competent person mandated by them. Do not modify or alter safety-relevant parts! The manufacturer cannot be held liable for any direct or indirect accident after a modification or use other than that provided in this manual. Who modifies or has modified, prepares or has prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market. Ensure that all the necessary precautions are taken in order to avoid hazards that can arise as the result of contact with blood or body fluids. The product shall be used only by a person trained and competent in its safe use. There must be ensured that medical conditions of the user do not affect his safety during normal use of the equipment or in case of emergency. In case of doubt, consult the doctor. •**Description of the product:** The basket stretcher is a device designed for the rescue and transport of patients. It can be used in all rescue operations in which it is impossible to intervene with the standard transport methods and where the patient needs protection from possible side impact. A sling system delivered with the stretcher. •**DESCRIPTION:** (Lanyard) The Stretcher lifting sling has been designed for rescue operation and is manufactured from raw material conforming to the requirements of fall arrest standards. It is specially designed for vertical transport of stretchers. The lanyards are connected to two length adjusting webbing legs to maintain the stretcher in the horizontal position. The four karabiners lock onto each corner of the stretcher to enable it to be lifted or lowered horizontally. Adjustment of lengths if required can be done by with the help of the buckle provided on the product. •**OPERATION and PRECAUTIONS:** Strap the patient properly to the stretcher before attempting any lift. Do a trial lift to check the level of the stretcher before attempting the final lift. Before and during use, we recommend that you make the necessary provisions for a possible rescue in complete safety. The length of the anchor lines and the types of equipment available in the workplace area are important factors. Risk assessment and rescue training should, therefore, equally be taken into account. The device must be checked each time before use so that any function faults and / or damage caused during transport and/or during storage are detected. Check the state of the lanyard before each use (no cuts in the wires) and make sure that the markings are readable. The lanyard should not be re-used if you have any doubts on the equipment state or after a fall. Excessive Rust. The stretcher lifting sling should not be re-used in the event of deformation, damage or doubt. It is forbidden to add or replace any lanyard component whatsoever. The assembly has a safe working load limit of:270 kg. Before and during use, we recommend that you make the necessary provisions for a possible rescue in complete safety. The length of the anchor lines and the types of equipment available in the workplace area are important factors. Risk assessment and rescue training should, therefore, equally be taken into account. Check the state of the lanyard before each use (no cuts in the wires) and make sure that the markings are readable. The lanyard should not be re-used if you have any doubts on the equipment state or after a fall. •This device can be used for helicopter transport. The stretcher can be used for HORIZONTAL and VERTICAL displacement. However, whatever the kind of displacement, it's important to maintain the stretcher in horizontal position. The stretcher is designed to be secured in the helicopter. •**DONNING AND/OR ADJUSTMENTS:** On receipt of the product: Remove packaging and layout all components so they are clearly visible. Check that all parts indicated on the shipping document are present. Check that the product is correctly assembled and that all the rivets are tight. Attach the safety belts by knotting them to the perimeter rope (fig. B, C) The device must be checked each time before use so that any function faults and / or damage caused during transport and/or during storage are detected. •**OPERATION:** Practice and training should be done with the stretcher without a patient so that the operators have familiarity with manoeuvring before regular use in emergency situations.. •**PATIENT LOADING:** Get as close as possible to the patient before loading. Load the trauma patient using the standard immobilising devices: spine board, vacuum mattress, etc.. •**FOOTREST:** Place the footrest at a sufficient height such that the planar surface is in contact with the patient's feet in order to prevent the longitudinal displacements. Insert the buckles on the end of the belts through the button holes/handles on the perimeter of the stretcher. •**LIFTING THE STRETCHER WITH A PATIENT ON BOARD:** Restrain the patient in the stretcher using the supplied belts. Check that the belts are not too tight. The operators must be positioned at the ends of the stretcher, one in head area the other by the feet. Using a correct lifting technique which will avoid excessive strain, each operator should take a strong grasp of the handholds at each end of the basket and lift. •**TRANSPORT ON INCLINED/SLOPING GROUND:** Correctly position the footrest so that the patient will not move.: If the patient has leg injuries, immobilise with the correct device and do not upright the stretcher.: Secure the stretcher with a static safety system before lifting or descending. The patient must always be assisted during manoeuvres. •**VERTICAL TRANSPORT:** Correctly position the footrest so that the patient will not move.. If

the patient has leg injuries, immobilise with the correct device and do not upright the stretcher. The patient must always be assisted during manoeuvres. **Usage limits:** •Read the instruction manual before use. Follow the approved EMS procedures for the restraint and transport of patients. Failure to do so could cause serious harm to the patient. This device is designed for use by at least 2 operators in good physical condition. Maximum load: 270 Kg The patient must always be immobilised on the stretcher. Do not use if the weight of the patient is not evenly distributed. Do not exceed the maximum loading weight and follow the maintenance guide indicated in this instruction manual. Do not leave the patient unassisted on the stretcher. The product should not be exposed or come into contact with sources of heat or inflammable agents. The operators must have the following minimum requirements: Have the physical strength to use the stretcher / Be able to tightly hold the stretcher with each hand / Have a generally good muscular condition.

**•CONTRAINDICATIONS AND SIDE EFFECTS:** The use of this device, if used as described in this manual, does not present any contra-indications or collateral effects. **•CONTROLS:** Before each use: Check that the perimeter rope is tight. Check that the safety belts are correctly attached and adapt for immobilising the patient. Make sure that the footrest is in position. If it has no defaults, it is therefore proper for the intended use. Contact us for further information; If you discover a defect, you should return the product to the retailer you purchased the product from or directly to DELTA PLUS GROUP. The device does not require periodical servicing. **Storage/Cleaning instructions:** Store and transport inside its original packaging and make sure it cannot be crushed, hit, dropped, or otherwise damaged by impact, and that it's not stored under heavy objects. Store in a cool dry place, away from light. Clean exposed parts with a sponge and delicate detergent then dry with a soft cloth or chammy leather. The mattress inside the stretcher can be cleaned using soap and water and disinfected using a normal disinfectant. Failure to carry out the cleaning procedure could result in cross contamination caused by residuals of secretions or body liquids. The correct and routine cleaning of the stretcher will guarantee it a longer life time. When the stretcher is put into service, a maintenance and check up routine should be arranged and a person responsible for its execution should be identified.

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TÜV RHEINLAND SHANGHAI LTD BUILDING 177 LANE 777 WEST GUANGZHONG ROAD 20007 SHANGHAI CHINE.